


ME 8SI



Technical data

Type		ME 8SI
Max. pumping speed 50/60 Hz according to ISO 21360	m ³ /h	7.2 / 7.8
Motor power	kW	0.2
Rated current at 100-120 V~ 50/60 Hz	A	4.0 / 4.2
200-230 V~ 50/60 Hz	A	2.0 / 2.1
Max. permitted range of supply voltage Attention: voltage selection switch		100-120 V~ +5/-10% 50-60 Hz 200-230 V~ +/-10% 50-60 Hz
Inlet		via protective filter
Dimensions L x W x H	mm	386 x 265 x 285
Weight	kg	16.2

Components	Wetted parts
Hose	PVC
Hose nozzle	PPS
O-ring	NBR
Valve block (adjustment device)	aluminium alloy
Seal ring manometer	copper
Hollow bolt, dispensing screw (adjustment device)	stainless steel